



MANUAL HOSPITAL BED BGY1J1 (BHBD-526)

MANUAL HOSPITAL BED BGY1J1

MANUAL BED FOUR CRANK 5 FUNCTION

Manual Hospital Bed is suitable for various operations conducted by medical units and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance



The main bed and auxiliary bed are connected with each other, auxiliary table is easy for dismantle. The auxiliary table with four wheels
The main bed adjusts various movements for the gears and auxiliary bed is ratched adjustment
Various Accessories are complete

SPECIFICATIONS

Model	BGY1J1
Old Model	BHBD-526
Table Length	1900±50 mm
Table width	600±20 mm
Height Adjustment	750±50 mm
Back board upwards	≥50°
Backboard downward	≥15°
Seat board downwards	≥8°
Seat board upwards	≥30°
Leg board upwards	≥30°
Leg board downwards	≥30°
Standard accessory	A pair of arm support, a pair of leg support, a pair of arm rods, a pair of foot pedal
Alt Name	Manual bed four crank 5 function

APPLICATIONS

Recuperating, Nursing



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com