



WHEELCHAIR BGY1K1 (BHBD-902)

WHEELCHAIR BGY1K1

MANUAL WHEELCHAIR

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.



It is a type of steady and simple Manual Wheelchair

SPECIFICATIONS

Model	BGY1K1
Old Model	BHBD-902
Wheelchair length	970 mm
Opened wheelchair width	660 mm
Wheelchair Height	920 mm
Closed wheelchair width	350 mm
Seat Width	430 mm
Front wheel	6" inch
Rear wheel	24" inch
Seat Height	470 mm
Seat depth	400 mm
Backrest Height	450 mm
Armrest Height	210 mm
Distance between armrest	425 mm
Height between footrests and ground	120 mm
Net Weight	12.94 Kg
Alt Name	Manual Wheelchair

APPLICATIONS

Patient Indoor-outdoor locomotion, Patient transfer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com