



## WHEELCHAIR BGY1O1 (BHBD-906)

# WHEELCHAIR BGY101

## MANUAL WHEELCHAIR

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.



It is a type of steady and simple Manual Wheelchair

## SPECIFICATIONS

Model	BGY101
Old Model	BHBD-906
Wheelchair length	1250 mm
Opened wheelchair width	690 mm
Wheelchair Height	1210 mm
Closed wheelchair width	365 mm
Seat Width	450 mm
Front wheel	8" inch
Rear wheel	24" inch
Seat Height	520 mm
Seat depth	455 mm
Backrest Height	450 mm
Armrest Height	240 mm
Distance between armrest	450 mm
Height between footrests and ground	100 mm
Net Weight	25.4 Kg
Alt Name	Manual Wheelchair

## APPLICATIONS

Patient Indoor-outdoor locomotion, Patient transfer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)