



WHEELCHAIR BGY1Q1 (BHBD-908)

WHEELCHAIR BGY1Q1

MANUAL WHEELCHAIR

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.



It is a type of steady and simple Manual Wheelchair

SPECIFICATIONS

Model	BGY1Q1
Old Model	BHBD-908
Wheelchair length	830 mm
Opened wheelchair width	580 mm
Wheelchair Height	780 mm
Closed wheelchair width	350 mm
Seat Width	360 mm
Front wheel	4" inch
Rear wheel	24" inch
Seat Height	470 mm
Seat depth	380 mm
Backrest Height	330 mm
Armrest Height	215 mm
Distance between armrest	415 mm
Height between footrests and ground	100 mm
Net Weight	13 Kg
Alt Name	Manual Wheelchair

APPLICATIONS

Patient Indoor-outdoor locomotion, Patient transfer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com