



LED OPERATION LIGHT BGZ-1802

LED OPERATION LIGHT BGZ-1802

DOUBLE ARM LED OPERATING LAMP

Challenging tasks require special solutions. The innovative LED operation light gives you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.



SPECIFICATIONS

Model	BGZ-1802
Illumination	40000~120000 LuX
Color Temperature (K)	3000~5000
Color Rendering Index (Ra)	≥98
Special Color Rendering Index (R9)	≥98
Special Color Rendering Index (R13/R15)	≥99
Diameter of Light Spot (mm)	120-350
Illumination Depth (mm)	≥1500
Lamp Head Power (VA)	180
LED Service Life (h)	≥600000 Hour(s)
Global Voltages	100-240V 50/60 HZ
Alt Name	Double Arm LED Operating Lamp

APPLICATIONS

Dental implant surgery, Medical unit, Surgery, Operation theatre.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com