



ELECTRIC HOSPITAL BED BGZ1C1

ELECTRIC HOSPITAL BED BGZ1C1

LUXURY FIVE FUNCTION BED WITH WEIGHING SYSTEM



SPECIFICATIONS

Model	BGZ1C1
Size	2200*1050*500/830mm
Back section adjustment	65°
Leg section adjustment	35°
Head tilt	≥12°
Foot tilt	≥15°
Height range	500-830mm
Safe Working Load	240 kg
One button to be chair function	Yes
0 reset function	Yes
Security automatic lock	Yes
Electric CPR	Yes
Alt Name	Luxury Five Function Bed With Weighing System

ITEMS INCLUDED





No	Name	Quantity
1	Bedhead board	2 pcs
2	ABS Guardrails	4 pcs
3	Central control	1 pc
4	Patient controller	2 pcs
5	JIECANG motors	5 pcs
6	Battery backup	1 pc
7	Angle indicator	4 pcs
8	IV hole	4 pcs
9	Uraine hook	4 pcs
10	Silient Universal castors	4 pcs
11	Bumpers	4 pcs

ACCESSORIES FOR PURCHASE

No	Name
1	Foot controller
2	Night light
3	Emergency stop button

4	Mattress
5	Overbed Food Table
6	Bedside cabinet
7	Color is optional
8	X-ray Bedboard
9	Weight scale function
10	Oxygen cylinder rack
11	Monkey bar

FEATURES



Weight
Weight System



90°

Angle Indicator



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com