



STRECHER BGZ1G1 (BHBD-814)

STRECHER BGZ1G1

HYDRAULIC TROLLEY CAN BE FILMED



SPECIFICATIONS

Model	BGZ1G1
Old Model	BHBD-814
Size	2140*810*630/930mm
Back angle	0°-90°
Leg angle	0°-30°
Load bearing	250kg
Regulation of tilt	0°-15°
Alt Name	Hydraulic Trolley Can Be Filmed

ITEMS INCLUDED

No	Name	Quantity
1	Crank handle	1 pc
2	X-ray back plate	1 pc
3	Footboard	3 pcs
4	Mattress	1 pc
5	Silient Universal castors	4 pcs

ACCESSORIES FOR PURCHASE

No	Name
1	Double opening guardrail
2	IV hole
3	Record table
4	Conductive wheel
5	Weighing function



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com