



## ELECTRIC HOSPITAL BED BGZ1I1

# ELECTRIC HOSPITAL BED BGZ1I1

## AUTOMATIC BED

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital beds allow height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.



With the nurse controller and CPR function  
Bed board adopts HPL board, x-ray photograph for whole body is available  
2 step adjustment at shank section

## SPECIFICATIONS

Model	BGZ1I1
Size	1950 x 910 x 440/800 mm
Height Adjustment	440-800 mm
Back-rest Adjustment	65°
Knee-rest Adjustment	35°
Trendelenburg	12°
Anti-Trendelenburg	12°
Bed Length Extension	2350 mm
Bed Panel Extension	Manual / 180 mm
Safe Working Load	240 kg
Motors	3 motors
Alt Name	Automatic Bed

## APPLICATIONS

Medical unit, ICU.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)