



MANUAL HOSPITAL BED BGZ1K1 (BHBD-524)

MANUAL HOSPITAL BED BGZ1K1

MANUAL HOSPITAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



Material: ABS

4 pcs of European style guardrail, upside down structure

Luxury central control castors, brake the bed with one stamp

With 4 bumpers

SPECIFICATIONS

Model	BGZ1K1
Old Model	BHBD-524
Size (L x W x H)	2280 x 1120 x 420/780 mm
Height Adjustment	420-780 mm
Back-rest Adjustment	65° (±10°)
Knee-rest Adjustment	35° (±10°)
Weight Loading Capacity	> 240 kg
Alt Name	Manual Hospital Bed

ITEMS INCLUDED

No	Name	Quantity
1	European Style Guardrail	4 pcs
2	ABS Head & Foot Board	2 pcs
3	5" Castors with Brake	4 pcs
4	Urine Hook	1 pc
5	Protective Bumpers	4 pcs
6	IV Pole	1 pc
7	Crank Quantity	3 pcs

APPLICATIONS

Recuperating, Nursing.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com