



MANUAL HOSPITAL BED BGZ1L1 (BHBD-525)

MANUAL HOSPITAL BED BGZ1L1

MANUAL HOSPITAL BED WITH 2 FUNCTION

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



SPECIFICATIONS

| | |
|----------------------|-------------------------------------|
| Model | BGZ1L1 |
| Old Model | BHBD-525 |
| Size (L x W x H) | 2200 x 1020 x 420/780 mm |
| Back-rest Adjustment | 65° (±10°) |
| Knee-rest Adjustment | 35° (±10°) |
| Height Range | 420-780 mm |
| Alt Name | Manual Hospital Bed With 2 Function |

ITEMS INCLUDED

| No | Name | Quantity |
|----|--------------------------|----------|
| 1 | European style guardrail | 4 pcs |
| 2 | ABS head and foot board | 2 pcs |
| 3 | Crank quantity | 2 pcs |
| 4 | 5" Castors with brake | 4 pcs |
| 5 | Urine hook | 1 pc |
| 6 | Protective bumpers | 4 pcs |
| 7 | IV pole | 1 pc |

APPLICATIONS

Recuperating, Nursing.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com