



STRECHER BGZ101 (BHBD-812)

STRECHER BGZ101

HOSPITAL TRANSFER STRETCHER TROLLEY

Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.



It is a type of manual transfer stretcher. Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down;

- Aluminum-Alloy protective rails, can be folded and hidden under the surface on no use;
- Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily;
- Rotate the crank and handle on the bedboard, the patient on the mattress can be moved parallelly to another bed or flat table.

SPECIFICATIONS

Model	BGZ101
Old Model	BHBD-812
Hi-Low range	920mm-610mm
Weight	89 kgs
Alt Name	Hospital Transfer Stretcher Trolley

APPLICATIONS

Patient transfer, Patient trolley.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com