



LED OPERATION LIGHT BGZ1Q1 (BOPL-406)

LED OPERATION LIGHT BGZ1Q1

SINGLE ARM LED OPERATING LAMP

Challenging tasks require special solutions. The innovative LED operation light gives you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.



SPECIFICATIONS

Model	BGZ1Q1
Old Model	BOPL-406
Technical Parameters	700 LED head
Illumination Intensity	40,000 - 160,000 Lux
Color Temperature	3000 - 5000 K
Color Rendering Index (Ra)	≥ 98
Special Color Rendering Index (R9)	≥ 98
Special Color Rendering Index (R13 / R15)	≥ 99
Diameter of Light Spot	120 - 350 mm
Illumination Depth	≥ 1500 mm
Lamp Head Power	180 VA
LED Service Life	≥ 60,000 hours
Global Voltage	100-240 V, 50/60 Hz
Alt Name	Single Arm LED Operating Lamp

APPLICATIONS

Medical unit, Surgery, Operation theatre.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com