

## STRECHER BGZ1W1 (BHBD-811)

## STRECHER BGZ1W1

### STAINLESS STEEL STRETCHER TROLLEY

Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

## SPECIFICATIONS

Model	BGZ1W1
Old Model	BHBD-811
Dimension	1950x600x780 mm
Alt Name	Stainless Steel Stretcher Trolley

## APPLICATIONS

Patient transfer, Patient trolley



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)