



## STRECHER BGZ1X1 (BHBD-813)

# STRECHER BGZ1X1

## PE DOCKING CART

Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.



Main frame of the flat cart is jet-molding carbon steel, solid and good looking  
Cart surface and guardrail is structured by one-step molding ABS material, solid and elegant

Adopting Separate design: back plate is adjusted by air spring, could be operated by only one hand, convenient to use

The crank of the cart (the screw is equipped with a two-way protection device) can adjust height of the cart surface from 600 to 830 mm

Equipped with international advanced central-control braking system, stable and reliable, the oriental castors are easy to operate

Meeting requirement from overall design concept, more humanized, ideal for transporting patients in hospital.

Product composition:

## SPECIFICATIONS

Model	BGZ1X1
Old Model	BHBD-813
Cart A	
- Size	1790x450 mm
- Height	570~850 mm
- Movements distance	280±20 mm (manual)
Cart B	
Cart C	
- Back adjustment	90°±5° (manual)
Cart ABC size	3575x770 mm
Alt Name	PE Docking Cart

## FEATURES

Main frame of the flat cart is jet-molding carbon steel, solid and good looking

Cart surface and guardrail is structured by one-step molding ABS material, solid and elegant

Adopting Separate design: back plate is adjusted by air spring, could be operated by only one hand, convenient to use

The crank of the cart (the screw is equipped with a two-way protection device) can adjust height of the cart surface from 600 to 830 mm

Equipped with international advanced central-control braking system, stable and reliable, the oriental castors are easy to operate

Meeting requirement from overall design concept, more humanized, ideal for transporting patients in hospital.

Product composition:

This product has A car, B and C three major components, including car A and car B, each has A screw and four control of castor, senior C car there are two big barriers and ABS engineering plastics of bed plate.

## APPLICATIONS

Patient transfer, Patient trolley



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)