



LIQUID DISPENSER BHA1C2 (BPIP-202)

LIQUID DISPENSER BHA1C2

ADJUSTABLE GLASS-INJECTION DISPENSER



Availability for autoclave of the whole device
A non-fragile, durable silicone tube made of special materials
Convenience of use with high efficiency, and low repetitive error, capable of greatly reducing time and costs in research and experiment
Solvent Resistance: general acid and alkali
Material: PP

SPECIFICATIONS

Model	BHA1C2
Old Model	BPIP-202
Bottle Volume	400 ml
Volume	0.1-1 ml
Test Volume	0.1 ml
Accuracy	±0.02
Qty	1
Unit	pcs
Package	10 pcs/carton
Gross weight (kg)	3.9
Dimensions (m ³)	0.52x0.28x0.29
Alt Name	Adjustable Glass-Injection Dispenser translucent glass



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com