



## LIQUID DISPENSER BHA1C3 (BPIP-203)

# LIQUID DISPENSER BHA1C3

## ADJUSTABLE GLASS-INJECTION DISPENSER



Availability for autoclave of the whole device  
A non-fragile, durable silicone tube made of special materials  
Convenience of use with high efficiency, and low repetitive error, capable of greatly reducing time and costs in research and experiment  
Solvent Resistance: general acid and alkali  
Material: PP

## SPECIFICATIONS

Model	BHA1C3
Old Model	BPIP-203
Bottle Volume	400 ml
Volume	0.5-5 ml
Test Volume	0.5 ml
Accuracy	±0.04
Qty	1
Unit	pcs
Package	10 pcs/carton
Gross weight (kg)	3.9
Dimensions (m <sup>3</sup> )	0.52x0.28x0.29
Alt Name	Adjustable Glass-Injection Dispenser translucent glass



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)