



## LIQUID DISPENSER BHA1C4 (BPIP-204)

# LIQUID DISPENSER BHA1C4

## ADJUSTABLE GLASS-INJECTION DISPENSER



Availability for autoclave of the whole device  
A non-fragile, durable silicone tube made of special materials  
Convenience of use with high efficiency, and low repetitive error, capable of greatly reducing time and costs in research and experiment  
Solvent Resistance: general acid and alkali  
Material: PP

## SPECIFICATIONS

|                              |  |
|------------------------------|--|
| Model                        | BHA1C4   |
| Old Model                    | BPIP-204   |
| Bottle Volume                | 500 ml   |
| Volume                       | 0.5-5 ml   |
| Test Volume                  | 0.5 ml   |
| Accuracy                     | ±0.04  |
| Qty                          | 1  |
| Unit                         | pcs  |
| Package                      | 10 pcs/carton                                    |
| Gross weight (kg)            | 4.2  |
| Dimensions (m <sup>3</sup> ) | 0.62x0.25x0.315                                  |
| Alt Name                     | Adjustable Glass-Injection Dispenser amber glass |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)