



LIQUID DISPENSER BHA1C7 (BPIP-207)

LIQUID DISPENSER BHA1C7

ADJUSTABLE GLASS-INJECTION DISPENSER



Availability for autoclave of the whole device
A non-fragile, durable silicone tube made of special materials
Convenience of use with high efficiency, and low repetitive error, capable of greatly reducing time and costs in research and experiment
Solvent Resistance: general acid and alkali
Material: PP

SPECIFICATIONS

Model	BHA1C7
Old Model	BPIP-207
Bottle Volume	1000 ml
Volume	1-10 ml
Test Volume	1 ml
Accuracy	±0.1
Qty	1
Unit	pcs
Package	10 pcs/carton
Gross weight (kg)	2.9
Dimensions (m ³)	0.66x0.315x0.335
Alt Name	Adjustable Glass-Injection Dispenser plastic



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com