



## CHILDREN MEDICAL BED BHBD-105

# CHILDREN MEDICAL BED BHBD-105

## ELECTRIC 3 FUNCTION CHILDREN MEDICAL BED

Children's medical beds built to the most current standards while also meeting the clinician's needs which allows easy care and disinfection. The bed is equipped with Trendelenburg and reverse Trendelenburg positions. Simple and modern design for safety and protection of the child.



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, pink, with patient name card holder, option: cartoon sticker

Full length powder coated side rail, drop type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking S.S foot pedal with casters dia. 5"

With 3 motors from Taiwan

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes.

## SPECIFICATIONS

Model	BHBD-105
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	480~720 mm (without mattress)
IV pole	Yes
Dimension	1980Lx980Wx480~720 mm
Alt Name	Electric 3 Function Children Medical Bed

## APPLICATIONS

Medical unit, Nursing



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)