



## WHEELCHAIR BHBD-1105

## WHEELCHAIR BHBD-1105

### ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.



Thickened aluminum alloy body, strong bearing capacity  
180° opening and closing, easy to realize the transfer of various scenes  
Electric motor, one-button control, easy lifting operation, convenient and quick displacement  
Soft and comfortable armrest, comfortable to use  
360-degree universal wheel, with locking mechanism, shifting is more secure  
PU leather seat cushion, waterproof, can be showered, hollow design, a toilet can be placed below, Easy to extract

## SPECIFICATIONS

Model	BHBD-1105
Material	Steel, PU
Maximum loading weight	150 kg
Lifting Range	Seat: 30 cm, from 40 cm to 70 cm.
Wheel	Two front wheels are with brake
Door width, chair can pass it	At least 58 cm
It suites for bed	Height of bed from 9 cm to 70 cm
Waterproof level	Yes, IP44, top grade waterproof
Control	Remote control
Handles	All foldable
Dimensions	55x65x73 cm
Rated power	96 W
Voltage	DC 24 V
Power supply	Battery, there is a button on battery to power it on/off
Alt Name	Electric lift patient transfer Wheelchair

## APPLICATIONS

Patient transfer, patient lift, toilet-bath chair, Patient Indoor-outdoor locomotion



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)