



DELIVERY TABLE BHBD-208

DELIVERY TABLE BHBD-208

DELIVERY TABLE 4002

The Delivery table is a flexible design suitable for conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Imported LINAK motor various actions are controlled by foot controller
The base seat and upright column covering are made of high quality 304 stainless steel which is anti corrosive and can be easily clean
External folding auxiliary table is for easy dismantling
Equipped with stainless steel sewage basin and other accessories
Foot treadle brake device which is easy for moving obstetric table

SPECIFICATIONS

Model	BHBD-208
Table Length	2000+50 mm
Table width	600+10 mm
Height Adjustment	(650-900)+50 mm
Trendelenburg/rev-trend	≥10°/≥12°
Backboard upward	≥50°
Leg board outward	≥90°
Working arealength/width	600x600 mm
Standard accessory	A pair of arm support, a pair of leg support, a pair of foot pedal, a pair of foot arm rods, a stainless steel sewage basin, a foot controller
Alt Name	Delivery table 4002

APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com