



## DELIVERY TABLE BHBD-210

# DELIVERY TABLE BHBD-210

## ELECTRIC GYNECOLOGICAL TABLE

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Adopt high density memory sponge and is once formed, with no seam  
Hydraulic lifting by foot pedal, other actions are controlled by pneumatic  
Imported Y shaped sealing ring and durable  
Stainless steel concealed auxiliary table  
The base is made up of high strength ABS  
Equipped with foot brake device which is easy for moving the table

## SPECIFICATIONS

Model	BHBD-210
Table Length	1300+50 mm
Table width	550+30 mm
Height Adjustment	(700 950)+50 mm
Trendelenburg/rev-trend	≥3°/≥10°
Dimension of auxiliary platform	(350x450)+20 mm
Backboard downward	≥70°
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin
Alt Name	Electric gynecological table

## APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)