



## DELIVERY TABLE BHBD-211

# DELIVERY TABLE BHBD-211

## ELECTRIC GYNECOLOGICAL TABLE

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Whole stainless steel hidden auxiliary table with flexible extension and retraction  
Equipped with hidden sewage basin  
With foot control, can adjust the movements easily  
Imported LINAK gear motor, without noise, have a long service life

## SPECIFICATIONS

Model	BHBD-211
Table Length	1450+50 mm
Table width	520+30 mm
Height Adjustment	(400 900)+50 mm
Trendelenburg/rev-trend	≥15°/≥25°
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin
Power Supply	220V 50Hz
Alt Name	Electric gynecological table

## APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)