



ELECTRIC HOSPITAL BED BHBD-401

ELECTRIC HOSPITAL BED BHBD-401

ELECTRIC BED

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.



Central control silent casters: 150 mm
5-pieces Aluminium alloy guard rail
Detachable PP head & foot board
Surface electrostatic powder coating
Bed framework made from cold-rolled steel plates

SPECIFICATIONS

| | |
|-------------------------|--------------|
| Model | BHBD-401 |
| Back-rest lifting angle | 0°~75° |
| Foot-rest lifting angle | 0°~45° |
| Trendelenburg | 0°~10° |
| Reverse trendelenburg | 0°~10° |
| Load Bearing | 250 kg |
| Bed Surface Height | 500~1000 mm |
| Patient surface | 900x1960 mm |
| Drainage bag holders | Yes |
| IV pole | Yes |
| Alt Name | Electric bed |

APPLICATIONS

Medical unit, ICU



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com