



ELECTRIC HOSPITAL BED BHBD-405

ELECTRIC HOSPITAL BED BHBD-405

MULTI-FUNCTIONAL ELECTRIC BED



SPECIFICATIONS

Model	BHBD-405
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2115 mm (L) * 1045 mm (W) * 477 mm - 745 mm
Bed board dimension	1940 mm (L) * 900 mm (W)
Backrest	0 - 65°
Knee gatch	0 - 40°
Trendelenburg	0 - 12°
Reverse Trendelenburg	0 - 12°
HI-LOW	745 mm - 477 mm
Net weight	139 KG
Safe load capacity	< 250 KG
Castors	360° swivel castors (Dia. 150 mm) with central brake
Battery backup	Battery backup for power outages
Standard accessory	I.V. Pole
Alt Name	Multi-functional Electric Bed



FEATURES





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com