



ELECTRIC HOSPITAL BED BHBD-407

ELECTRIC HOSPITAL BED BHBD-407

MULTI-FUNCTIONAL ELECTRIC BED



Advanced durable & long-lasting heat-sensing touch panels.

The inner protective railing touch panels allow patients to access bed control and adjust as needed for their level of comfort. The outer protective railing touch panels allow caregivers to access bed control.

A 10.2" touch screen on the foot board can control all bed functions or lock them for safety purposes.

The postures can be saved and retrieved. This function helps nurses to enable specific posture quickly and easily.

Weighing scale and weight monitor functions. The weight of patients can be recorded and displayed as a graph. If patient's weight change exceeds the pre-set alarm line, the system will alarm with beep and signal.

The PP head & foot boards are integrally blow-molded, easy to detach and clean.

SPECIFICATIONS

Model	BHBD-407
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2280 mm (L) * 1050 mm (W) * 500 mm - 750 mm
Bed board dimension	1940 mm (L) * 900 mm (W)
Backrest	0 - 65°
Knee gatch	0 - 40°
Trend/Reverse Trend	0 - 12°
HI-LOW	750 mm - 500 mm
Net weight	158 KG
Safe load capacity	< 250 KG
Castors	360° swivel castors (Dia. 150 mm) with electric brake
Alt Name	Multi-functional Electric Bed



FEATURES

Advanced durable & long-lasting heat-sensing touch panels.

The inner protective railing touch panels allow patients to access bed control and adjust as needed for their level of comfort. The outer protective railing touch panels allow caregivers to access bed control.

A 10.2" touch screen on the foot board can control all bed functions or lock them for safety purposes.

The postures can be saved and retrieved. This function helps nurses to enable specific posture quickly and easily.

Weighing scale and weight monitor functions. The weight of patients can be recorded and displayed as a graph. If patient's weight change exceeds the pre-set alarm line, the system will alarm with beep and signal.

The PP head & foot boards are integrally blow-molded, easy to detach and clean.

4 damping PP protective railings with built-in angle indicators, controlled by gas springs. When protective railings fall and retract under the bed board, they can achieve seamless connection with the transfer stretcher.

Extendable abdomen & knee sections on bed board.

One touch electric brake.

IV pole sockets (one on each corner).

Battery backup for power outages.

Standard accessory: IV pole.

FUNCTIONS:

Backrest, Knee gatch, Height adjustment, Trend./Reverse Trend., Backrest and knee simultaneously moving, Reset, Weighing scale, Off-bed alarming, Electric brake.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com