



ELECTRIC HOSPITAL BED BHBD-408

ELECTRIC HOSPITAL BED BHBD-408

ELECTRIC 2 FUNCTION MEDICAL BED

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital beds allow height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.



ABS head & foot boards can be removable. 8kg/pair, color option: brown, blue, with patient name card holder
2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock
Bed sheet top made from perforated steel, epoxy powder coated and baked finish
Bed framework made from mild steel ERW rectangular tubes
With casters dia. 5" with individual brake
With 2 motors from Taiwan
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes
Optional Accessories: Shoe shelf, mosquito net rack, dinner table.

SPECIFICATIONS

Model	BHBD-408
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2120Lx970Wx500H mm
Alt Name	Electric 2 Function Medical Bed

APPLICATIONS

Medical unit, ICU



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com