



## ELECTRIC HOSPITAL BED BHBD-411

# ELECTRIC HOSPITAL BED BHBD-411

## ELECTRIC 2 FUNCTION MEDICAL BED

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital beds allow height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder  
ABS side rail, 4-piece, tuck away type with safe lock. with angle display  
Bed sheet top made from perforated steel, epoxy powder coated and baked finish  
Bed framework made from mild steel ERW rectangular tubes  
ABS foot pedal with casters dia. 5"  
With 2 Linak motors from Denmark. with battery  
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes.

## SPECIFICATIONS

Model	BHBD-411
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2260Lx1060Wx500H mm
Alt Name	Electric 2 Function Medical Bed

## APPLICATIONS

Medical unit, ICU.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)