



## ELECTRIC HOSPITAL BED BHBD-414

# ELECTRIC HOSPITAL BED BHBD-414

## ELECTRIC 3 FUNCTION MEDICAL BED

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital beds allow height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.



ABS head & foot boards can be removable. With bumper wheels, 12 kg/pair, color option: brown, blue, green with patient name card holder  
ABS side rail, 4pcs, tuck away type with safe lock. With angle display  
Bed sheet top made from perforated steel, epoxy powder coated and baked finish  
Bed framework made from mild steel ERW rectangular tubes  
Central locking both side foot pedal with casters dia. 5"  
With 3 linak motors from Denmark. With a battery.  
Backrest & knee-rest up/down at same time  
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes  
Optional Accessories: CPR/electric CPR.

## SPECIFICATIONS

Model	BHBD-414
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	430-720 mm (without mattress)
IV pole	Yes
Dimension	2260Lx1060Wx430-720 mm
Alt Name	Electric 3 function medical Bed

## APPLICATIONS

Medical unit, ICU.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)