



ELECTRIC HOSPITAL BED BHBD-419

ELECTRIC HOSPITAL BED BHBD-419

FIVE FUNCTION ELECTRIC HOSPITAL BED



SPECIFICATIONS

Model	BHBD-419
Size	2200*1100*420/780mm
Back-rest Adjustment	65°(±10°)
Knee-rest Adjustment	35°(±10°)
Trendelenburg	0-12°
Reverse Trendelenburg	0-12°
Alt Name	Five Function Electric Hospital Bed


ITEMS INCLUDED

No	Name	Quantity
1	JIECHANG Motors	4 pcs
2	Central Control	1 pcs
3	Hand Controller	1 pcs
4	Silent Universal Castors	4 pcs
5	ABS Guardrails	4 pcs

ACCESSORIES FOR PURCHASE


No	Name
1	Battery backup
2	X-ray bedboard
3	Mattress
4	Overbed Food Table
5	Bedside Cabinet

FEATURES



35°(±10°)

Leg-rest




420-780mm

Height
Adjustment




0-12°

Trendelenburg
angle



0-12°

Re-Trendelenburg
angle



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com