



MANUAL HOSPITAL BED BHBD-516

MANUAL HOSPITAL BED BHBD-516

MANUAL 3 FUNCTION MEDICAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



ABS head & foot boards can be removable. with bumper, 7kg/pair, color blue, with patient name card holder
2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock
Bed sheet top made from perforated steel, epoxy powder coated and baked finish
Bed framework made from mild steel ERW rectangular tubes
With casters dia. 5" with individual brake
Manual 3 functions adjusted by hidden type hand crank
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes
Optional Accessories: Shoe shelf, mosquito net rack, dinner table.

SPECIFICATIONS

Model	BHBD-516
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	430~720 mm
IV pole	Yes
Dimension	2120Lx980Wx430~720 mm
Alt Name	Manual 3 Function Medical Bed

APPLICATIONS

Recuperating, Nursing.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com