



## MANUAL HOSPITAL BED BHBD-519

# MANUAL HOSPITAL BED BHBD-519

## MANUAL 2 FUNCTION MEDICAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder  
ABS side rail, 4-piece, tuck away type with safe lock. with angle display  
Bed sheet top made from perforated steel, epoxy powder coated and baked finish  
Bed framework made from mild steel ERW rectangular tubes.  
Central locking both side foot pedal with casters dia. 5"  
Manual 2 functions adjusted by hidden type hand crank  
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes.

## SPECIFICATIONS

Model	BHBD-519
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2220Lx1060Wx500 mm
Alt Name	Manual 2 Function Medical Bed

## APPLICATIONS

Recuperating, Nursing.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)