



MANUAL HOSPITAL BED BHBD-521

MANUAL HOSPITAL BED BHBD-521

FLAT HOSPITAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



SPECIFICATIONS

Model	BHBD-521
Size	2150*900*500mm
Material	ABS plastic with punching
Alt Name	Flat Hospital Bed

APPLICATIONS

Recuperating, Nursing



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com