



MANUAL HOSPITAL BED BHBD-523

MANUAL HOSPITAL BED BHBD-523

CRANKS HOSPITAL BED



SPECIFICATIONS

Model	BHBD-523
Size	2150*960*500mm
Back section adjustment	$\geq 75^\circ$
Loading capacity	$\geq 240\text{kg}$
Uraine hook	Yes
Bumpers	Yes
Alt Name	Cranks Hospital Bed

ITEMS INCLUDED

No	Name	Quantity
1	Stainless steel handle	1 pc
2	Aluminum alloy Guardrails	2 pcs
3	Universal 125 castors with braking	4 pcs
4	Bedhead board	2 pcs
5	IV Pole	1 pc

ACCESSORIES FOR PURCHASE

No	Name
1	Night light
2	Mattress
3	Overbed Food Table
4	Bedside cabinet
5	Color is optional

FEATURES





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com