



## MANUAL HOSPITAL BED BHBD-525

# MANUAL HOSPITAL BED BHBD-525

## MANUAL HOSPITAL BED WITH 2 FUNCTION

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



## SPECIFICATIONS

Model	BHBD-525
Size (L x W x H)	2200 x 1020 x 420/780 mm
Back-rest Adjustment	65° (±10°)
Knee-rest Adjustment	35° (±10°)
Height Range	420-780 mm
Alt Name	Manual Hospital Bed With 2 Function

## ITEMS INCLUDED

No	Name	Quantity
1	European style guardrail	4 pcs
2	ABS head and foot board	2 pcs
3	Crank quantity	2 pcs
4	5" Castors with brake	4 pcs
5	Urine hook	1 pc
6	Protective bumpers	4 pcs
7	IV pole	1 pc

## APPLICATIONS

Recuperating, Nursing.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)