



## STRECHER BHBD-810

# STRECHER BHBD-810

## PATIENT STRETCHER

Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.



## SPECIFICATIONS

Model	BHBD-810
Frame Material	Thickened 202 stainless steel
Mattress Material	Leather (integrated leather pad)
Mattress Size	1750-500 mm, thickness 4 cm
Product Size	1900 x 600 x 750 mm
Wheel Features	Universal anti-skid wheels (two with brakes, two without brakes)
Wheel Size	5 inches
Product Color	Black
Alt Name	Patient Stretcher

## APPLICATIONS

Patient transfer, Patient trolley.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)