

STRECHER BHBD-811

STRECHER BHBD-811

STAINLESS STEEL STRETCHER TROLLEY

Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

SPECIFICATIONS

Model	BHBD-811
Dimension	1950x600x780 mm
Alt Name	Stainless Steel Stretcher Trolley

APPLICATIONS

Patient transfer, Patient trolley



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com