



## STRECHER BHBD-812

# STRECHER BHBD-812

## HOSPITAL TRANSFER STRETCHER TROLLEY

Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.



It is a type of manual transfer stretcher. Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down;

- Aluminum-Alloy protective rails, can be folded and hidden under the surface on no use;
- Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily;
- Rotate the crank and handle on the bedboard, the patient on the mattress can be moved parallelly to another bed or flat table.

## SPECIFICATIONS

Model	BHBD-812
Hi-Low range	920mm-610mm
Weight	89 kgs
Alt Name	Hospital Transfer Stretcher Trolley

## APPLICATIONS

Patient transfer, Patient trolley.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)