



STRECHER BHBD-818

STRECHER BHBD-818

HYDRAULIC TRANSFER STRETCHER



1. A transfer stretcher with double hydraulic systems. Functions operated by treading the foot pedals at both sides.
2. Clear arrow symbols to instruct the user to operate.
3. Central lockable 360° swivel castors (Dia. 150mm). Retractable 5th wheel provides effortless directional movement and cornering.

SPECIFICATIONS

Model	BHBD-818
Overall dimension	1906*743 mm (L*W)
Height range (bed board to ground)	813 - 590 mm
Bed board dimension	1906*610 mm (L*W)
Backrest	0 - 85°
Trend/Reverse Trend	0 - 12°
Protective railings	Damping PP protective railings, lifting controlled by gas spring, seamless connection with stretcher/operation table
Net Weight	82 KG
Standard accessories	Mattress, I.V. Pole
Alt Name	Hydraulic Transfer Stretcher

FEATURES

1. A transfer stretcher with double hydraulic systems. Functions operated by treading the foot pedals at both sides.
2. Clear arrow symbols to instruct the user to operate.
3. Central lockable 360° swivel castors (Dia. 150mm). Retractable 5th wheel provides effortless directional movement and cornering.





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com