



STRECHER BHBD-819

STRECHER BHBD-819

HYDRAULIC TRANSFER STRETCHER



A transfer stretcher with double hydraulic systems. Six functions are adjusted by hydraulic system, gas spring and crank.
2 sets of push handles at the head and foot sides help users move the stretcher easily.
Central lockable 360° swivel castors (Dia. 150 mm). Retractable 5th wheel provides effortless directional movement and cornering.

SPECIFICATIONS

Model	BHBD-819
Overall dimension	2050*745 mm (L*W)
Height range (bed board to ground)	860 - 630 mm
Bed board dimension	1932*630 mm (L*W)
Backrest	0 - 85°
Knee gatch	0 - 50°
Trend/Reverse Trend	0 - 13°
Side rails	Foldable stainless steel side rails
Bed board	Radio-translucent bed board enables X-ray imaging on stretcher
Net Weight	105 KG
Standard accessory	I.V. Pole
Optional feature	Dia. 200 mm castors optional (Height range changes to: 910 - 680 mm)
Alt Name	Hydraulic Transfer Stretcher

FEATURES

A transfer stretcher with double hydraulic systems. Six functions are adjusted by hydraulic system, gas spring and crank.
2 sets of push handles at the head and foot sides help users move the stretcher easily.
Central lockable 360° swivel castors (Dia. 150 mm). Retractable 5th wheel provides effortless directional movement and cornering.





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com