



TISSUE EMBEDDING SYSTEM BHM1BD1 (BHTP-604)

TISSUE EMBEDDING SYSTEM BHM1BD1

PARAFFIN DISPENSER



- 1.The full automatic procedure control&new heating element gives feature of quick heating&energy saving
- 2.The temperature survey integration block from American DALLAS Corporation has been used to make it precise in temperature survey and credible in performance
- 3.It possesses the function of memorization to keep the set temperature automatically after startup
- 4.The LED shows heating state

SPECIFICATIONS

Model	BHM1BD1
Old Model	BHTP-604
Range of temperature	Ambient to 70 °C
Precision	±1 %
Voltage & Power	AC220V 50Hz / AC110V 60Hz, Power: 1200VA
Capacity	10000 ml
Dimension	350x450x510 mm
Net Weight	10 kg
Package Info	
--Packing	Wooden Case
--Gross Weight	15 kg
--Package Dimensions	510x450x530 mm
Alt Name	Paraffin Dispenser



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com