



TISSUE EMBEDDING SYSTEM BHM1BE1 (BHTP-606)

TISSUE EMBEDDING SYSTEM BHM1BE1

TISSUE WATER BATH



1. Model Number :HS1125
2. Short Desc :Tissue Water Bath
3. Heat conduct fast, durable.
4. Self memory of temperature.
5. Display presetting temperature & current temperature.
6. Corrosion-resistance of the bath.

SPECIFICATIONS

Model	BHM1BE1
Old Model	BHTP-606
Range of setting temperature	Ambient to 90 °C (constant adjusting)
Temperature precise error	±2 °C
Allowance ambient temperature	0-40 °C
Working voltage	220V / 50Hz
Power	200 W
Dimensions (WxDxH)	340x340x120 mm
Net Weight	4.7 kg
Package Info	
--Packing	Wooden Case
--Gross Weight	7 kg
--Package Dimensions	450x420x200 mm
Alt Name	Tissue Water Bath



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com