



## HOTPLATE MAGNETIC STIRRER BHMS-109

# HOTPLATE MAGNETIC STIRRER BHMS-109

## ECONOMICAL MAGNETIC STIRRER

This is magnetic stirrer customized for liquid mixing, suitable for chemical analysis, biology reagent mixing, regular laboratory experiments and liquid handling.



Maintenance free brushless DC motor  
Rare earth magnets give powerful stirring  
Stepless speed regulation 300-2000 rpm  
Ultra silence and vibration-free design  
PET top provides excellent chemical resistance  
Shaft can be fitted with a pH sensor or other sensors for different purposes

## SPECIFICATIONS

Model	BHMS-109
Work plate size	Ø120mm
Work plate material	PET
Motor type	Brushless DC motor
Stirring position	1
Max. stirring quantity (H2O)	1.5L
Max. magnetic bar (length)	40mm
Speed range	300-2000rpm
Speed setting	scale
Protection class	IP21
Motor output power	1.2W
Power	7.2W
Voltage, Frequency	100-240V, 50/60Hz or DC 12V
Dimension [WxDxH]	150x145x50mm
Weight	0.3kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	Economical Magnetic Stirrer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)