



## MICROPLATE READER BHN1B1 (BMRW-102)

# MICROPLATE READER BHN1B1

## MICROPLATE READER

Microplate Reader is a laboratory instrument that is used to measure chemical, biological or physical reactions, properties and analytes within the well of a microplate. Flexibility and reliability makes this product an unique choice.



Bichromatic, 8-channel Optical system.  
Large storage capacity up to 200 test protocols and 100,000 results.  
LED light source with lifespan over 10,000 hours.  
3 shaking modes with variable time.  
Windows operation system, self-check, malfunction alarm.  
7 inch touch TFT LCD screen, test results simultaneously displayed in screen.

## SPECIFICATIONS

|                         |   |
|-------------------------|---|
| Model                   | BHN1B1                                    |
| Old Model               | BMRW-102                                  |
| Repeatability           | <1%                                       |
| Stability               | <0.005 Abs/10 min                         |
| Linearity               | R <sup>2</sup> ≥0.995                     |
| Display                 | 7 inch TFT LCD, 800x480 Pixels            |
| Interface               | RS~232, 4 USB ports, SD Card              |
| Printer                 | Built-in printer, external printer        |
| Power Supply            | AC 100~240V, 50/60 Hz                     |
| Dimension               | 400Lx310Wx169H mm                         |
| Weight                  | 6 kg                                      |
| Readout Range           | 0.000~3.500 Abs                           |
| Plate Type              | 96-well plate                             |
| Light Source            | Monochromatic LED                         |
| Reading Speed           | 5s with single wavelength for whole plate |
| Input                   | Mouse, membrane keyboard, touch screen    |
| Humidity                | 15%~85%                                   |
| Temperature             | 10°C~35°C                                 |
| Measure System          | 8~channel optical system                  |
| Band Width and Accuracy | 8 nm ± 2 nm                               |
| Shaking                 | 3 shaking modes with variable time        |
| Wavelength              | 405, 450, 492, 630nm, two more optional   |
| Alt Name                | Microplate Reader                         |

## APPLICATIONS

Drug Discovery, Research, Laboratory, Medical, Assay Development, Life Science.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)