



## TISSUE EMBEDDING SYSTEM BHTP-607

# TISSUE EMBEDDING SYSTEM BHTP-607

TISSUE BAKING PROCESSOR



1. Heat conduct fast, durable.
2. Self memory of temperature.
3. Display presetting temperature & current temperature.
4. Corrosion-resistance of the bath.

## SPECIFICATIONS

Model	BHTP-607
Range of setting temperature	Ambient to 90 °C (constant adjusting)
Temperature precise error	±2 °C
Allowance ambient temperature	0-40 °C
Working voltage	220V / 50Hz
Power	250 W
Dimensions (WxDxH)	340x340x120 mm
Net Weight	5.3 kg
Package Info	
--Packing	Wooden Case
--Gross Weight	7 kg
--Package Dimensions	450x400x200 mm
Alt Name	Tissue Baking Processor



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)