



TISSUE EMBEDDING SYSTEM BHTP-609

TISSUE EMBEDDING SYSTEM BHTP-609

TISSUE FLOATING WATER BATH



The round tissue bath has black non-stick coated interior to reduce glare and improve visibility. The top ledge allows space for preheating or for evaporating away any remaining water. A durable fiberglass insulated body prevents heat loss. Glass lid and detachable power cord for easy cleaning.

SPECIFICATIONS

Model	BHTP-609
Capacity	2.5 L
Temperature range	RT to 75 °C
Outside Size	270 mm (D) x 210 mm (H)
Interior Size	210 mm (D) x 65 mm (H)
Power consumption	400 W
Power supply	AC 110V/220V ±10%, 50/60Hz
Net Weight	3 kg
Package Info	
--Packing	Wooden Case
--Gross Weight	4 kg
--Package Dimensions	340x350x160 mm
Alt Name	Tissue Floating Water Bath



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com