



## BULLET ICE MAKER BIBU-101

# BULLET ICE MAKER BIBU-101

VERTICAL ICE MAKER



Luxurious appearance and compact design;  
Fully automatically computer controlled, convenient for use;  
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;  
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;  
Cylindrical bullet ice with elegant appearance.  
Ice thickness can be adjusted, with a wider scope of application.  
Bullet Ice

## SPECIFICATIONS

Model	BIBU-101
Ice making capacity	20kg/24h
Power/voltage	220w
Ice storage capacity	6kg
Ice shape	Bullet shape
N.W./kgs	28kg
G.W./kgs	31kg
Product dimension	330x470x605mm
Packing Dimension	400x550x680
20GP	168
40GP	348
40HQ	464
Alt Name	Vertical Ice Maker

## FEATURES

Luxurious appearance and compact design;  
Fully automatically computer controlled, convenient for use;  
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;  
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;  
Cylindrical bullet ice with elegant appearance.  
Ice thickness can be adjusted, with a wider scope of application.



Bullet Ice



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)