



BULLET ICE MAKER BIBU-104

BULLET ICE MAKER BIBU-104

VERTICAL ICE MAKER




Luxurious appearance and compact design;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Cylindrical bullet ice with elegant appearance.
Ice thickness can be adjusted, with a wider scope of application.
Bullet Ice

SPECIFICATIONS

| | |
|----------------------|--------------------|
| Model | BIBU-104 |
| Ice making capacity | 50kg/24h |
| Power/voltage | 380w |
| Ice storage capacity | 14kg |
| Ice shape | Bullet shape |
| N.W./kgs | 35kg |
| G.W./kgs | 38kg |
| Product dimension | 400x510x825mm |
| Packing Dimension | 480x590x890mm |
| 20GP | 98 |
| 40GP | 196 |
| 40HQ | 290 |
| Alt Name | Vertical Ice Maker |

FEATURES

Luxurious appearance and compact design;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Cylindrical bullet ice with elegant appearance.
Ice thickness can be adjusted, with a wider scope of application.



Bullet Ice



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com