



## BULLET ICE MAKER BIFL-212

# BULLET ICE MAKER BIFL-212

VERTICAL ICE MAKER



Luxurious appearance and compact design;  
Fully automatically computer controlled, convenient for use;  
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;  
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;  
Cylindrical bullet ice with elegant appearance.  
Ice thickness can be adjusted, with a wider scope of application.  
Bullet Ice

## SPECIFICATIONS

Model	BIFL-212
Ice making capacity	120kg/24h
Power/voltage	760w
Ice storage capacity	50kg
Ice shape	Bullet shape
N.W./kgs	70kg
G.W./kgs	75kg
Product dimension	680x610x1050mm
Packing Dimension	750x680x1100mm
20GP	48
40GP	96
40HQ	96
Alt Name	Vertical Ice Maker

## FEATURES

Luxurious appearance and compact design;  
Fully automatically computer controlled, convenient for use;  
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;  
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;  
Cylindrical bullet ice with elegant appearance.  
Ice thickness can be adjusted, with a wider scope of application.



Bullet Ice



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)