



CUBE ICE MAKER BIFL-216

CUBE ICE MAKER BIFL-216

SPLIT-TYPE ICE MAKER



Split-type structure, large ice bin, large ice storage capacity and strong applicability;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Square ice cubes with beautiful shapes.
Square

SPECIFICATIONS

Model	BIFL-216
Ice Making Capacity	1000kg/24h
Ice Storage Capacity	500kg
Ice Shape	Square
Power/Voltage (W/V)	4500w/220v
Cooling Type	Water cooling
Product Dimension (mm)	1250x950x2000
Machine Appearance	Split machine
Alt Name	Split-Type Ice Maker

FEATURES

Split-type structure, large ice bin, large ice storage capacity and strong applicability;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Square ice cubes with beautiful shapes.



Square



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com